

pCAGGS

PVT1020 2ug

Lot no 1901093

Store at -20°C

Expiry date

pCAGGS Information

Promoter: CAG promoter

Replicator: SV40 ori, ori

Terminator: beta -globin poly (A) signal

Plasmid classification: lactation serial plasmid; lactation expression plasmid; pCAG series plasmid.

Plasmid size: 4801bp

Prokaryotic resistance: ampicillin Amp (100 u g/ml)

Cloned strains of Escherichia coli, DH5 A and other Escherichia coli

Culture conditions: 37 centigrade, aerobic, LB

Expression host: mammalian cells such as 293T

Culture conditions: 37 °C, 5%CO₂

Induction mode: no induction, instantaneous expression

5'sequencing primers: pcaggs-F (GTTCGGCTTCTGGCGTGT)


3'sequencing primers: pcaggs-R (TATGTCCTCCGAGTGAGAG)

pCAGGS Description

pCAGGS plasmid can be used to express gene efficiently under the control of chicken β -actin, rabbit β -globulin heterozygous promoter (CAG), and human CMV-IE enhancer in various mammalian cells. The CAG promoter sequence is part of the chicken β -actin promoter, the first exon and the first intron (it seems to have strong enhanced subtype activity. " It is linked to the rabbit β -globin fragment, the 3' part of the second intron includes the branching point needed for normal splicing reaction and the 5' portion of the third exon.

This plasmid is useful for highly efficient expression of genes under the control of the AG promoter and the human CMV-IE enhancer in various mammalian cells. The AG promoter sequence consists of the chicken β -actin promoter, the first exon and part of the first intron (that seems to have a strong enhancer-like activity) linked to a rabbit β -globin fragment, consisting of a 3' part of the second intron (inclusive a branch point which is required for normal splicing reactions) and a 5' part of the third exon. When cloning a fragment downstream from the lac promoter it may be advisable to use lacIq strains in order to prevent fortuitous expression of a possibly noxious polypeptide.

pCAGGS Map



source 1..4801
/organism="synthetic DNA construct"
/mol_type="other DNA"

polyA_site 158..213
/note="beta-globin poly(A) signal"
/note="rabbit beta-globin polyadenylation signal"

primer_bind complement(574..590)
/note="M13 rev"
/note="common sequencing primer, one of multiple similar variants"

protein_bind 598..614
/bound_moiety="lac repressor encoded by lacI"
/note="lac operator"
/note="The lac repressor binds to the lac operator to inhibit transcription in E. coli. This inhibition can be relieved by adding lactose or isopropyl-beta-D-thiogalactopyranoside (IPTG)."

promoter complement(622..652)
/note="lac promoter"
/note="promoter for the E. coli lac operon"

protein_bind 666..687
/bound_moiety="E. coli catabolite activator protein"
/note="CAP binding site"

/note="CAP binding activates transcription in the presence
of cAMP."
promoter 746..941
/note="SV40 promoter"
/note="SV40 early promoter"
rep_origin 792..927
/note="SV40 ori"
/note="SV40 origin of replication"
polyA_signal 947..1081
/note="SV40 poly(A) signal"
/note="SV40 polyadenylation signal"
rep_origin complement(1319..1907)
/direction=LEFT
/note="ori"
/note="high-copy-number ColE1/pMB1/pBR322/pUC origin of
replication"
CDS complement(2078..2938)
/codon_start=1
/gene="bla"
/product="beta-lactamase"
/note="AmpR"
/note="confers resistance to ampicillin, carbenicillin, and
related antibiotics"

Manual

RVGYI /translation="MSIQHFRVALIPFFAAAFCLPVFAHPETLVKVKDAEDQLGA
QNDLVEYS ELDLNSGKILESFRPEERFPMMSTFKVLLCGAVLSRIDAGQEQLGRRRIHYS
DHVTRLDRW PVTEKHLTDGMTVRELCSAAITMSDNTAANLLLTIGGPKELTAFLHNMG
VAGPLLRSA EPELNEAIPNDERDTTMPVAMATTLRKLLTGELLTLASRQQLIDWMEADK
RQIAEIGAS LPAGWFIADKSGAGERGSRGIIAALGPDGKPSRIVVIYTTGSQATMDERN

LIKHW"

promoter complement(2939..3043)

/gene="bla"

/note="AmpR promoter"

enhancer 3074..3453

/note="CMV enhancer"

/note="human cytomegalovirus immediate early enhancer"

promoter 3455..3732

/note="chicken beta-actin promoter"

intron 3733..4750

/note="chimeric intron"

/note="chimera between introns from chicken beta-actin and
rabbit beta-globin"

ORIGIN

1 ttctcgagg aattcactcc tcaggtgcag gctgcctatc agaaggtggt ggctggtgtg

61 gccaatgccc tggctcacia ataccactga gatcttttc cctctgcaa aaattatggg
121 gacatcatga agccccttga gcatctgact tctggctaataaaggaatt tattttcatt
181 gcaatagtgt gttggaattt tttgtgtctc tactcggaa ggacatatgg gagggcaaat
241 catttaaac atcagaatga gtatttggtt tagagtttg caacatatgc ccatatgctg
301 gctgcatga acaaagggtg gctataaaga ggcatcagt atatgaaaca gcccctgct
361 gtccattcct tattccatag aaaagccttg acttgaggtt agatttttt tatattttg
421 tttgtgttat tttttcttt aacatcccta aaatttcct tacatgtttt actagccaga
481 ttttctcc tctcctgact actcccagtc atagctgtcc ctcttctt atggagatcc
541 ctgacctgc agccaagct tggcgtaatc atggcatag ctgttctctg tgtgaaattg
601 ttatccgctc acaattcac acaacatag agccggaagc ataaagtga aagcctgggt
661 gcctaatgag tgagtaact cacattaatt gcgttgcgct cactgcccgc tttccagtcg
721 ggaaacctgt cgtgccagcg gatccgcatc tcaattagtc agcaaccata gtcccggccc
781 taactccgcc catcccggcc ctaactccgc ccagtccgc ccatttccg ccccatggct
841 gactaatttt tttatttat gcagaggccg aggccgctc ggctctgag ctattccaga
901 agtagtgagg aggtttttt ggaggcctag gctttgcaa aaagctaact tgtttattgc
961 agcttataat ggttacaat aaagcaatag catcacaat tcacaaata aagcattttt
1021 ttcactgcat tctagtgtg gttgtccaa actcatcaat gtatctatc atgtctggat
1081 ccgctgcatt aatgaatcgg ccaacgcgcg gggagaggcg gtttgcgtat tggcgctct
1141 tccgcttct cgtcactga ctcgctgcgc tcggtcgtc ggctgcggcg agcggtatca
1201 gctcactcaa aggcggtaat acggttatcc acagaatcag ggataacgca ggaagaaca
1261 tgtgagcaaa aggcagcaa aaggccagga accgtaaaaa ggccgcttg ctggcgtttt
1321 tccataggct ccgccccct gacgagcatc acaaaaatcg acgctcaagt cagagggtgc
1381 gaaaccggac aggactataa agataccagg cgtttcccc tgaagctcc ctcgtcgcgt

1441 ctctgttcc gaccctgccg cttaccggat acctgtccgc ctttctcct tcggaagcg
1501 tggcgcttc tcatagctca cgctgtagg atctcagttc ggtgtaggtc gttcgctcca
1561 agctgggctg tgtgcacgaa cccccgttc agcccgaccg ctgcgcctta tccgtaact
1621 atcgtcttga gtccaaccg gtaagacacg acttatcgcc actggcagca gccactggta
1681 acaggattag cagagcgagg tatgtaggcg gtgctacaga gttcttgaag tggtagccta
1741 actacggcta cactagaaga acagtattg gtatctgcgc tctgctgaag ccagttacct
1801 tcggaaaag agttgtagc tcttgatccg gcaaacaac caccgctgtt agcgggtggt
1861 ttttgtttg caagcagcag attacgcgca gaaaaaagg atctcaagaa gatcctttga
1921 tctttctac ggggtctgac gctcagtga acgaaaactc acgttaaggg atttggta
1981 tgagattatc aaaaaggatc ttcacctaga tccttttaa ttaaaatga agtttaaat
2041 caatctaaag tatatatgag taaacttgg ctgacagtta ccaatgctta atcagtgagg
2101 cacctatctc agcgatctgt ctatttcgtt catccatagt tgctgactc cccgtcgtg
2161 agataactac gatacgggag ggcttaccat ctggccccag tgctgcaatg ataccgcgag
2221 acccagctc accggctcca gatttatcag caataacca gccagccgga agggccgagc
2281 gcagaagtgg tctgcaact ttatccgct ccatccagtc tattaattgt tgccgggaag
2341 ctgagtaag tagttcgcca gtaaatagtt tgcgcaacgt tgttgccatt gctacaggca
2401 tcgtggtgc acgctcgtc tttggatgg cttcattcag ctccggtcc caacgatcaa
2461 ggcgagttac atgatcccc atgtgtgca aaaaagcgg tagctcttc ggtcctccga
2521 tcgtgtcag aagtaagttg gccgcagtgt taccatcat ggttatggca gcactgcata
2581 attctctac tgcatgcca tccgtaagat gctttctgt gactggtgag tactcaacca
2641 agtcattctg agaatagtt atgcggcgac cgagttgctc ttgccggcg tcaatcggg
2701 ataataccgc gccacatagc agaacttaa aagtgtcat cattggaaa cgttcttcgg
2761 ggcgaaaact ctcaaggatc ttaccgctgt tgagatccag ttcgatgtaa cccactcgtg

2821 cacccaactg atcttcagca tcttttactt tcaccagcgt tctgggtga gcaaaaacag
2881 gaaggcaaaa tgccgcaaaa aaggaataa gggcgacacg gaaatgtga atactcatac
2941 tcttcctttt tcaatattat tgaagcattt atcagggtta ttgtctcatg agcggataca
3001 tatttgaatg tatttagaaa aataaacaaa taggggttcc ggcacattt ccccgaaaag
3061 tgccacctgg gtcgacattg attattgact agtattaat agtaatcaat tacggggta
3121 ttagttcata gccatataat ggagttccgc gttacataac ttacggtaaa tggcccgcct
3181 ggctgaccgc ccaacgacc cgcgccattg acgtcaataa tgacgtatgt tcccatagta
3241 acgccaatag ggactttcca ttgacgtcaa tgggtggagt atttacggta aactgccac
3301 ttggcagtac atcaagtga tcatatgcca agtacgccc ctattgacgt caatgacgt
3361 aaatggccc cctggcatta tgcccagtac atgacctat gggacttcc tacttggcag
3421 tacatctacg tattagtcac cgctattacc atggctgagg tgagccccac gttctgctc
3481 actctcccca tctcccccc ctccccacc ccaattttgt atttattat ttttaatta
3541 ttttgtcag cgatgggggc gggggggggg gggggggccc cccagggcg ggcggggcg
3601 ggcgaggggc ggggcggggc gaggcggaaa ggtgcggcg cagccaatca gagcggcg
3661 ctccgaaagt ttcctttat ggcgagggcg cggcggcggc ggcctataa aaagcgaagc
3721 gcgcgggcg cgggagtcgt tgcgctgc cttcccccg tgccccgctc cgccgccgc
3781 tcgcccgc cgccccgct ctgactgacc gcgttactcc cacaggtgag cggcgggac
3841 ggccttctc ctccgggctg taattagcgc ttggttaat gacggctgt tcttttctg
3901 tggctgcgtg aaagcctga ggggctccg gagggcctt tgtcggggg gagcggctc
3961 ggggtgcgt gcgtgtgtg gtgcgtggg agcggcgct gcggctccg gctcccggc
4021 ggctgtgagc gctcggggc cggcggggg cttgtgcgc tccgactgt gcgagggg
4081 agcggggcg gggcggtgc cccggtgc gggggggct gcgagggaa caaaggctg
4141 gtgcggggtg tgtgcgtgg ggggtgagca ggggtgtgg gcgctcgt cgggctgca

Manual

4201 cccccctgc accccctcc ccgagttgct gagcacggcc cggcttcggg tgcggggctc
4261 cgtacggggc gtggcgctgg gctcgccgtg ccgggcgggg ggtggcggca ggtgggggtg
4321 ccgggcgggg cggggccgcc tcgggcccgg gagggctcgg gggaaggggc gcggcggccc
4381 ccggagcgcc ggcggctgtc gaggcgctgc gagccgcagc cattgccttt tatgtaatc
4441 gtgcgagagg ggcagggac ttctttgtc ccaaactgtg gcggagccga aatctgggag
4501 gcgccccgc acccctcta gcgggcgctg ggcgaagcgg tgcggcggc gcaggaagga
4561 aatgggcggg gagggccttc gtgcgtgcc gcgccccgt ccccttctcc ctctccagcc
4621 tcggggctgt ccgcgggggg acggctgcct tcggggggga cggggcaggg cggggttcgg
4681 ctctggcgt gtgaccggcg gctctagagc ctctgtaac catgtcatg ccttctctt
4741 tttctacag ctctgggca acgtgctggt tattgtgctg tctcatcatt ttggcaaga
4801 a

//

Caution:

1. This product is FOR RESEARCH USE ONLY!