

Human Hepatocellular carcinoma-associated protein TD26, TD26 Recombinant Protein

Packing:1mg

Purity:>95%

Search name

Human betatrophin Recombinant Protein, Human Hepatocellular carcinoma-associated protein TD26, TD26 Recombinant Protein, Human TD26 Recombinant Protein, Human C19orf80 Recombinant Protein, Human ANGPTL8 Recombinant Protein, Human PRO1185 Recombinant Protein, Human PVPA599 Recombinant Protein, Human RIFL Recombinant Protein, Human hepatocellular carcinoma-associated gene TD26 Recombinant Protein, Human lipasin Recombinant Protein, Human chromosome 19 open reading frame 80 Recombinant Protein

Accession: NP_061157.3 (22aa-198aa)

Gene name: C19orf80

Gene Accession: NM_018687.6

Host: E coli

MW: 36.47 KD

Storage buffer: Tris-based Buffer(PH8.0) with 10% glycerol

Storage buffer: Tris-based Buffer(PH8.0) with 10% glycerol

12% SDS-PAGE:



