

**Human Major facilitator superfamily domain- containing protein 2A, MFSD2A Protein
For research use only**

Description: recombinant protein with His-tag (partial)

Packing: 1mg

Purity: >90%

Search name

Human MFSD2A PROTEIN ,Human MFSD2 PROTEIN ,Human major facilitator superfamily domain containing 2 PROTEIN ,Human major facilitator superfamily domain-containing protein 2A PROTEIN ,Human major facilitator superfamily domain containing 2A PROTEIN ,Human MFSD2A PROTEIN ,Human HMFN0656 PROTEIN ,Human PP9177 PROTEIN ,Human UNQ300/PRO341 PROTEIN ,Human MFSD2 PROTEIN, Human Sodium-dependent lysophosphatidylcholine symporter 1 PROTEIN

Protein Accession: NP_116182.2 (453aa-530aa)

Sequence:

agyqtrgcsqpervkftlnmlvtmapivlillglllflkmpideerrrqnkkaqlrdeasssgcsetdstelasil

Gene name: MFSD2A

Gene Accession: NM_032793.4

Host: E coli

MW: 25.58 KD

Applications: AP, Array, ELISA, WB-Re

Storage buffer: Tris-based Buffer(PH8.0) with 10% glycerol

Storage Instruction: Store at -80°C. (avoid repeated freezing and thawing.)

12% SDS-PAGE:

